


Product Datasheet

Human Cystatin M/Cystatin E/CST6 ELISA Kit (orb526882)

Description	Human Cystatin M/Cystatin E/CST6 ELISA Kit
Reactivity	Human
Range	46.9 pg/ml - 3,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: R29-M149
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA).
Uniprot ID	Q15828
Sensitivity	10 pg/ml
Dilution Range	Dilution ratio of 1:5, concentration in serum and plasma is around 3 ng/ml.
Expiration Date	Please enquire.



Human
Cystatin
M/CST6
PicoKine
ELISA Kit...