


Product Datasheet

Rat Serpin D1 ELISA Kit (orb526885)

Description	Rat Serpin D1 ELISA Kit
Reactivity	Rat
Range	0.78 ng/ml - 50 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: E24-S479
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (EDTA, citrate) and urine.
Uniprot ID	Q64268
Sensitivity	0.1 ng/ml
Dilution Range	Dilution ratio of 1:200, concentration in serum and plasma is around 1.2 ug/ml.
Expiration Date	Please enquire.



Rat Serpin
D1 PicoKine
ELISA
Kitstandard...