

Product Datasheet

BAT3/BAG6 Antibody (orb527016)

Description

BAT3/BAG6 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, FC, ICC, IF, IHC, WB

Applications

Immunogen

E. coli-derived human BAT3/BAG6 recombinant protein (D15-A88).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Tested Species: In-house tested species with positive results.
By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

150 kDa

Uniprot ID


[P46379](#)

Dilution Range


Western blot, 0.1-0.5µg/ml, Human, Rat
Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, Rat, By Heat Immunohistochemistry (Frozen Section), 0.5-1µg/ml, Human
Immunocytochemistry/Immunofluorescence, 2µg/ml, Human
Flow Cytometry, 1-3µg/1x10⁶ cells, Human

Expiration Date


12 months from date of receipt.



Flow Cytometry analysis of SH-SY5Y cells...



WB analysis of BAT3/BAG6 using anti-BAT3...



IF analysis of BAT3/BAG6 using anti-BAT3...