



Product Datasheet

HECTD3 Antibody (orb527025)


Description	HECTD3 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-Fr, WB
Immunogen	A synthetic peptide corresponding to a sequence of human HECTD3 (HYASAKVCEEKLRYYAAYNCVAIDTDMSPWEE).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	97 kDa
Uniprot ID	Q5T447
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A549 cells US...



WB analysis of HECTD3 using anti-HECTD3 ...



IHC analysis of HECTD3 using anti-HECTD3...