



Product Datasheet

RMI2 Antibody (orb527037)


Description	RMI2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC, IHC-Fr, WB
Immunogen	A synthetic peptide corresponding to a sequence of human RMI2 (KMTDLSDNPIHESMWEELEVEDLHRNIP).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19 kDa
Uniprot ID	Q96E14
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunohistochemistry (Frozen Section), 0.5-1µg/ml Immunocytochemistry, 0.5-1µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Flow
Cytometry
analysis of
U2OS cells
US...



WB analysis of
RMI2 using
anti-RMI2
anti...



IHC analysis of
RMI2 using
anti-RMI2
ant...