


Product Datasheet

ZP3 Antibody (orb527042)

Description	ZP3 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	A synthetic peptide corresponding to a sequence of human ZP3 (LRLMEENWNAEKRSPFHLGDAHLQAEIHT).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	47 kDa
Uniprot ID	P21754
Dilution Range	Western blot, 0.1-0.5µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of ZP3 using anti-ZP3 antibo...