



Product Datasheet

Cystatin C/CST3 Antibody (monoclonal, 4H8) (orb527056)


Description	Cystatin C/CST3 Antibody (monoclonal, 4H8)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E. coli-derived human Cystatin C recombinant protein (Position: K31-A146).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	15 kDa
Uniprot ID	P01034
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 2µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells ELISA (Cap), 1-5µg/ml
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A549 cells us...



WB analysis of Cystatin C/CST3 using ant...



IF analysis of CST3 using anti-CST3 anti...