



Product Datasheet

DHFR Antibody (monoclonal, 3C8) (orb527060)


Description	DHFR Antibody (monoclonal, 3C8)
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	E.coli-derived human DHFR recombinant protein (Position: V2-D187). Human DHFR shares 90% amino acid (aa) sequence identity with both mouse and rat DHFR.
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	22 kDa
Uniprot ID	P00374
Dilution Range	Western blot, 0.1-0.5 µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1 µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of DHFR using anti-DHFR anti...



IHC analysis of DHFR using anti-DHFR ant...



IHC analysis of DHFR using anti-DHFR ant...