



Product Datasheet

GM130 GOLGA2 Antibody (monoclonal, 4G3) (orb527062)


Description	GM130 GOLGA2 Antibody (monoclonal, 4G3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human GM130 recombinant protein (Position: E796-E913).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	140 kDa
Uniprot ID	Q08379
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 5µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A431 cells us...



WB analysis of GM130 using anti-GM130 an...



IF analysis of GM130 using anti-GM130 an...