



Product Datasheet

STUB1 Antibody (monoclonal, 1318) (orb527067)


Description	STUB1 Antibody (monoclonal, 1318)
Species/Host	Mouse
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC, WB
Immunogen	E.coli-derived human STUB1 recombinant protein (Position: R51-Y303).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
MW	35 kDa
Uniprot ID	Q9UNE7
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Flow Cytometry, 1-3µg/1x10 ⁶ cells
Expiration Date	12 months from date of receipt.



Flow Cytometry analysis of A549 cells using anti-STUB1 antibody.



WB analysis using anti-STUB1 antibody. Lane 1: Molecular Weight Marker; Lane 2: Blank; Lane 3: IgG; Lane 4: 1318; Lane 5: 1318 + BSA.



IHC analysis of STUB1 using anti-STUB1 antibody.