

Product Datasheet

Recombinant Human FGF2 Protein (orb527074)

Description	0
Reactivity	Human
Form/Appearance	Lyophilized
Storage	Lyophilized recombinant mouse Fibroblast Growth Factor basic (FGF-basic/FGF-2) remains stable up to 6 months at -80°C from date of receipt. Upon reconstitution, rhFGF2 remains stable up to 2 weeks at 4°C or up to 3 months at -20°C.
Note	For research use only
Protein Sequence	PALPEDGGSGAFPPGHFKDPKRLYCKNGGFFLRHPDGRVDGVREKSDPHIKLQ LQAEERGVSIGVVCANRYLAMKEDGRLLASKCVTDECFFFERLESNNYNTYRS RKYTSWYVALKRTGQYKLGSKTGPGQKAILFLPMSAKS
Source	E. Coli-derived P143-S288
MW	16.5KD
Uniprot ID	P09038
Expiration Date	6 months from date of receipt.