

Product Datasheet

ATG12 antibody (orb527085)

Description	Mouse monoclonal antibody to ATG12
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	ATG12
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b, Igk
Clonality	Monoclonal
MW	15113
Uniprot ID	O94817
NCBI	NP_001264712.1 , NP_004698.3
Dilution Range	WB: 1:4000
Expiration Date	12 months from date of receipt.