

Product Datasheet

CAB39 antibody (orb527086)

Description	Mouse monoclonal antibody to CAB39
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	CAB39
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1,K
Clonality	Monoclonal
MW	39869
Uniprot ID	Q9Y376
NCBI	NP_001124322.1 , NP_001124321.1 , NP_057373.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.