

Product Datasheet

ENO2 antibody (orb527088)

Description	Mouse monoclonal antibody to ENO2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Recombinant Protein
Target	ENO2
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG2b,k
Clonality	Monoclonal
MW	47269
Uniprot ID	P09104
NCBI	NP_001966.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.