

Product Datasheet

UBE2C antibody (orb527107)

Description	Mouse monoclonal antibody to UBE2C
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthesized Peptide
Target	UBE2C
Form/Appearance	Purified monoclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein G column, followed by dialysis against PBS.
Note	For research use only
Isotype	IgG1, κ
Clonality	Monoclonal
MW	19652
Uniprot ID	000762
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.