

Product Datasheet

CXCR3 antibody (orb527119)

Description	Rabbit polyclonal antibody to CXCR3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	CXCR3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40660
Uniprot ID	P49682
NCBI	NP_001136269.1, NP_001495.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FC: 1:10-50
Expiration Date	12 months from date of receipt.