

Product Datasheet

SETD1B antibody (orb527146)

Description	Rabbit polyclonal antibody to SETD1B
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	SETD1B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	212803
Uniprot ID	Q9UPS6
NCBI	NP_055863.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.