

Product Datasheet

BRD9 antibody (orb527176)

Description	Rabbit polyclonal antibody to BRD9
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	BRD9
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	67000
Uniprot ID	Q9H8M2
NCBI	NP_001009877.2 , NP_076413.3
Dilution Range	WB: 1:1000- 1:2000
Expiration Date	12 months from date of receipt.