



Product Datasheet

NEUROG3 antibody (orb527185)


Description	Rabbit polyclonal antibody to NEUROG3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	NEUROG3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23077
Uniprot ID	Q9Y4Z2
NCBI	NP_066279.2
Dilution Range	IHC-P: 1:25, WB: 1:2000, IF: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of N...



Western blot: All lanes: Anti-Neurogenin3 ...



Western blot: Anti-Neurogenin3 Antibody (...