

Product Datasheet

ERBB4 antibody (orb527186)

Description	Rabbit polyclonal antibody to ERBB4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	ERBB4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	146808
Uniprot ID	Q15303
NCBI	NP_001036064.1 , NP_005226.1
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.