

Product Datasheet

MORC2 antibody (orb527189)

Description	Rabbit polyclonal antibody to MORC2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	MORC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	117823
Uniprot ID	Q9Y6X9
Dilution Range	WB: 1:2000, IF: 1:25
Expiration Date	12 months from date of receipt.