

Product Datasheet

MAP4K1 antibody (orb527202)

Description	Rabbit polyclonal antibody to MAP4K1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Target	MAP4K1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	91296
Uniprot ID	Q92918
Dilution Range	WB: 1:1000-2000
Expiration Date	12 months from date of receipt.