



Product Datasheet

IDUA antibody (orb527236)


Description	Rabbit polyclonal antibody to IDUA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	IDUA
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	72670
Uniprot ID	P35475
NCBI	NP_000194.2
Dilution Range	FC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of lysates from A5...



Western blot: Anti-IDUA Antibody (Center)...



Western blot: All lanes: Anti-IDUA Antibod...