



Product Datasheet

PKM antibody (orb527246)


Description	Rabbit polyclonal antibody to PKM
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	PKM
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	57937
Uniprot ID	P14618
NCBI	NP_001193728.1 , NP_001193727.1 , NP_001193726.1 , NP_872270.1 , NP_872271.1 , NP_002645.3 , NP_001193725.1
Dilution Range	WB: 1:2000, IF: 1:200, IHC-P: 1:50- 1:100
Expiration Date	12 months from date of receipt.



Fluorescent confocal image of MCF7 cells...



Western blot: All lanes: Anti-PKM2 Antibod...



Western blot analysis of hPKM2-N491 in A...