

Product Datasheet

E2F1 antibody (orb527255)

Description	Rabbit polyclonal antibody to E2F1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	E2F1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46920
Uniprot ID	Q01094
NCBI	NP_005216.1
Dilution Range	WB: 1:500
Expiration Date	12 months from date of receipt.