

Product Datasheet

KIR3DL1 antibody (Biotin) (orb52726)

Description	KIR3DL1 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 3DL1 protein (81-336AA)
Target	KIR3DL1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P43629
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.