

Product Datasheet

PTGER3 antibody (orb527263)

Description	Rabbit polyclonal antibody to PTGER3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	PTGER3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43310
Uniprot ID	P43115
NCBI	NP_001119516.1 , NP_942012.1 , NP_942009.1 , NP_942011.1 , NP_942007.1 , NP_942010.1 , NP_942008.1
Dilution Range	WB: 1:1000, FC: 1:10-50
Expiration Date	12 months from date of receipt.