

Product Datasheet

TNFSF11 antibody (orb527302)

Description	Goat polyclonal antibody to TNFSF11
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF
Target	TNFSF11 / OPGL
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Notes: Immunofluorescence: Strong expression of the protein seen in the membranes of MCF7 and NIH3T3 cells. Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of MCF7 cells. Recommended concentration
Protein Sequence	DLAKRSKLEAQP
Clonality	Polyclonal
MW	35.5; 27.7
Entrez	8600
Dilution Range	ELISA: 1:16000, IHC-P: 10ug/ml
Expiration Date	12 months from date of receipt.