

Product Datasheet

Dehydroepiandrosterone-BSA conjugate [Out of stock] (orb527315)

Description	Dehydroepiandrosterone-BSA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C.Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Carrier protein: Bovine serum albumin (BSA)
CAS Number	53-43-0
Expiration Date	12 months from date of receipt.