

Product Datasheet

Dehydroepiandrosterone-OVA conjugate (orb527318)

Description	Dehydroepiandrosterone-OVA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L,pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C.Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Other Small Molecule Antigen
Source	Ovalbumin (OVA)
CAS Number	53-43-0
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.