

Product Datasheet

Gibberellin-BSA conjugate (orb527328)

Description	Gibberellin-BSA conjugate
Tested Applications	ELISA, IHC
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C.Long-term store at -23°C. Avoid repeated freeze.
Note	For research use only
Application notes	Other Small Molecule Antigen
Source	Bovine serum albumin (BSA)
CAS Number	77-06-5
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.