

Product Datasheet

Salbutamol-BSA conjugate (orb527342)

Description	Salbutamol-BSA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C, Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Food Safety Small Molecule Antigen
Source	Bovine serum albumin (BSA)
CAS Number	18559-94-9
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.