

Product Datasheet

The metabolite of Nitrofurantoin-BSA conjugate (orb527360)

Description	The metabolite of Nitrofurantoin-BSA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L,pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C.Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Food Safety Small Molecule Antigen
Source	Bovine serum albumin (BSA)
CAS Number	67-20-9
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.