

Product Datasheet

Retinoic acid- HSA conjugate (orb527386)

Description	Retinoic acid- HSA conjugate
Tested Applications	ELISA, ICA
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C.Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Other Small Molecule Antigen
Source	Human serum albumin (HSA)
CAS Number	302-79-4
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.