

Product Datasheet

S-Adenosyl-L-homocysteine -BSA conjugate (orb527388)

Description	S-Adenosyl-L-homocysteine -BSA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C. Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Other Small Molecule Antigen
Source	Bovine serum albumin (BSA)
CAS Number	979-92-0
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.