

Product Datasheet

Chloroamphenicol-BSA conjugate (orb527400)

Description	Chloroamphenicol-BSA conjugate
Tested Applications	ELISA, GICA, ICA
Preservatives	0.01mol/L, pH 7.4 PBS
Form/Appearance	Dry powder
Storage	Short-term store at 4°C, Long-term store at -20°C. Avoid repeated freeze.
Note	For research use only
Application notes	Food Safety Small Molecule Antigen
Source	Bovine serum albumin (BSA)
CAS Number	56-75-7
Hazard Information	For research use only
Expiration Date	12 months from date of receipt.