

Product Datasheet

Goat IgG Fab Fragment protein (orb527417)

Description	Goat IgG Fab Fragment protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Goat IgG Fab fragment
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Goat serum IgG digested with papain
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.