

Product Datasheet

Streptolysin O Protein (orb527425)

Description	Streptolysin O Protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Streptolysin O protein
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from Streptococcus hemolyticus
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.