

Product Datasheet

Human Myoglobin Protein (orb527430)

Description	Human Myoglobin Protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	human Myoglobin protein
Preservatives	20 mM Tris-HCl, 0.5 M NaCl, pH 8.0
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from Human heart tissue
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.