

Product Datasheet

Chicken Yolk Immunoglobulin protein (orb527456)

Description	Chicken Yolk Immunoglobulin protein
Tested Applications	ELISA, SDS-PAGE, WB
Target	Chicken Yolk Immunoglobulin
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Aliquot and store at -20°C or -80°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Purity	> 95%(SDS-PAGE)
Source	Purified from chicken egg yolk
Hazard Information	For research use only
Expiration Date	6 months from date of receipt.