



Product Datasheet

MBP antibody (orb52747)


Description	MBP antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Myelin basic protein (1-197AA)
Target	MBP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P02686
Dilution Range	WB:1:500-5000, IHC-P:1:100-500, IF:1:50-500
Expiration Date	12 months from date of receipt.



Western blot analysis
of Rat brain(lane ...



Immunohistochemical
staining of human
br...



Immunofluorescence
analysis of Hela cell...