

Product Datasheet

PCCA antibody (FITC) (orb52829)

Description	PCCA antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Propionyl-CoA carboxylase alpha chain, mitochondrial protein (585-720AA)
Target	PCCA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05165
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.