

Product Datasheet

PNLIPRP2 antibody (HRP) (orb52868)

Description	PNLIPRP2 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Pancreatic lipase-related protein 2 protein (305-460AA)
Target	PNLIPRP2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P54317
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.