

Product Datasheet

GABBR1 (phospho-Ser923) antibody (orb5288)

Description	GABBR1 (phospho-Ser923) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GABBR1 around the phosphorylation site of Ser923 (R(p-S)RR)
Target	phospho-GABBR1 (Ser923)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBS5
Expiration Date	12 months from date of receipt.