



Product Datasheet

PRKAR2A antibody (orb52907)


Description	PRKAR2A antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human cAMP-dependent protein kinase type II-alpha regulatory subunit protein (43-205AA)
Target	PRKAR2A
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P13861
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa(lane 1), M...



Immunohistochemical staining of human ly...



Immunofluorescence analysis of HepG2 cel...