



Product Datasheet

PSMC5 antibody (orb52943)

Description	PSMC5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human 26S proteasome regulatory subunit 8 protein (253-324AA)
Target	PSMC5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P62195
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...



Immunofluorescence analysis of PC3 cells...