

Product Datasheet

ROCK2 antibody (FITC) (orb52969)

Description	ROCK2 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human
Conjugation	FITC
Immunogen	Recombinant Human Rho-associated protein kinase 2 protein (1109-1388AA)
Target	ROCK2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075116
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.