

Product Datasheet

GATA1 (phospho-ser142) antibody (orb5302)

Description	GATA1 (phospho-ser142) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human GATA1 around the phosphorylation site of ser142 (RL(p-S)PD)
Target	phospho-GATA1 (ser142)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P15976
Expiration Date	12 months from date of receipt.