



Product Datasheet

WISP3 antibody (orb53095)


Description	WISP3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human WNT1-inducible-signaling pathway protein 3 protein (201-372AA)
Target	WISP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O95389
Dilution Range	WB:1:500-5000, IHC-P:1:20-200
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 whole cel...



Immunohistochemical staining of human li...



Immunohistochemical staining of human ki...