

Product Datasheet

acta1b antibody (Biotin) (orb53154)

Description	acta1b antibody (Biotin)
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Zebrafish Actin protein (1-168AA)
Target	acta1b
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6XNL8
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.